Simpleware Medical Instructions for Software Use

Version V-2024.06

Last amended: July 2024 (revision number 2.7).

Table of Contents

ENGLISH	2
ČEŠTINA	8
DANSK	14
ESPAÑOL	20
FRANÇAIS	26
ITALIANO	32
MAGYAR	38
SUOMI	44
SVENSKA	50

©2024 Synopsys, Inc. All rights reserved. Synopsys is a trademark of Synopsys, Inc. in the United States and other countries. A list of Synopsys trademarks is available at <http://www.synopsys.com/copyright.html>. All other names mentioned herein are trademarks or registered trademarks of their respective owners.

Simpleware Medical Instructions for Software Use

Version V-2024.06

Last amended: July 2024 (revision number 2.7).

This document is an electronic copy of the basic 'Instructions for Use' of Simpleware Medical Software. Detailed instructions can be found in the 'Reference Guide'.

1 General Description of Simpleware Medical

Simpleware™ Medical is standalone software package for converting 3D image data into models that can be exported in different file formats for applications as per the intended use.

2 Intended Use

Simpleware Medical is intended for use as a software interface and image segmentation system for the transfer of medical imaging information to an output file. The software is intended to inform clinical management by improving confidence in clinical decision-making and pre-surgical planning through use of patient-specific 3D models (virtual 3D models and 3D printed anatomical models) created from medical imaging. Outputs from the software must be reviewed for clinical suitability by expert clinicians working on specific patient cases.

2.1 Clinical Benefits

Simpleware Medical provides potential indirect benefits to intended users:

- Improved understanding of individual patient anatomies through 3D model visualization and segmentation as an improvement on 2D visualization
- Improved planning of surgeries through virtual 3D models and 3D printed anatomical models
- Improved ability to select safest/optimal treatment (implant choice/device sizing and customization to specific patients) through accurate measurements of anatomy
- Improved future clinical procedures through analysis of past procedures by comparing pre-surgical plans to post-clinical imaging
- Improved clinical training and patient communication through patient-specific 3D models

Simpleware Medical provides potential direct clinical benefits for patients through improved treatment planning:

- Lower risk of complications and improved recovery times
- Improved functional and aesthetic surgical results, enhancing patient's psychological and physical recovery
- Reduced likelihood of a need for revision surgery

2.2 Indications for Use

Simpleware Medical is intended for all medical applications where medical imaging is used to aid in diagnosis or treatment planning. There are no limitations for disease type, tissue, and/or part of the body, and patient population. Please note that the medical device is indicated for these uses except where contraindicated.

2.3 Patient Population

Simpleware Medical has no intended patient population. Outputs (virtual 3D models and 3D printed anatomical models) from the software can be used for informing clinical management for these criteria:

- Weight: No limitations in software
- Height: No limitations in software
- Health: No limitations in software
- Patient state: Patient is not user

2.4 Contraindications

Simpleware Medical is not intended for grading/sizing of anatomical defects to determine acceptable treatment plans. The output of the software is not a surgical plan itself, but helps with informing clinical management decisions – the software is not intended to drive clinical management decisions, or to directly treat clinical situations or patient conditions. In addition, Simpleware Medical is:

- Not intended for direct diagnosis or monitoring of vital physiological processes where the nature of variations of those parameters is such that it could result in immediate danger to the patient, for instance variations in cardiac performance, respiration, activity of the central nervous system, or diagnosis in clinical situations where the patient is in immediate danger.
- Not intended to provide inputs to directly treat or diagnose or drive clinical management for non-serious to critical scenarios (including death or irreversible deterioration of a person's state of health, serious deterioration of a person's state of health or a surgical intervention).
- Not intended to control or monitor the performance of active devices, or direct influence on the performance of these devices:
 - Active therapeutic class IIb devices
 - Active devices intended to emit ionizing radiation for therapeutic purposes
 - Active implantable devices
 - Active devices intended to emit ionizing radiation for diagnostic or therapeutic radiology, including interventional radiology devices
 - Devices used for direct diagnosis or monitoring of vital physiological processes where the nature of variations of those parameters is such that it could result in immediate danger to the patient, for instance variations in cardiac performance, respiration, activity of the central nervous system, or diagnosis in clinical situations where the patient is in immediate danger.
- Not intended for monitoring of physiological processes, including monitoring of vital physiological parameters, where the nature of variations of those parameters is such that it could result in immediate danger to the patient.
- Not intended to drive or influence the use of a device used to take decisions with diagnosis or therapeutic purposes where such decisions have an impact that may cause death or irreversible deterioration of a person's state of health, or a serious deterioration of a person's state of health or a surgical intervention.

2.5 Intended Users

The intended users of Simpleware Medical are biomedical engineers with a professional level of English. Specific training on the device is not required before use, but tutorials and a reference guide are available within the software if needed.

2.6 Intended Use Environment

Simpleware Medical is intended for use in a general office environment and workstation, assuming that intended user criteria are met. Please see **Supported Operating Systems** and **Recommended System Requirements** in this IFU for more details.

3 Technical Description

3.1 Warnings and Recommendations

The intended users of Simpleware Medical are biomedical engineers with a professional level of English. Specific training on the device is not required before use, but tutorials and a reference guide are available within the software if needed. English is the language used in the Simpleware Medical software interface.

The output must be verified by the responsible clinician.

No acceptance and performance testing, recurrent testing and maintenance is required for use. Software updates and patches will be provided by the manufacturer, where needed.

The DICOM Standard (version 3.0; 2022) is followed. Compliance details are available on the Simpleware website (<https://www.synopsys.com/simpleware/resources/regulatory.html>). We strongly recommend running through tests of the PACS export functionality on a PACS test system before using it in a clinical setting.

Simpleware Medical has the ability to process, store or discard information contained in medical image files such as DICOM files during the import process of these files. The import process can involve different data transfer methods including USB, CD/DVD, disk drive or server-based storage systems such as PACS. When such files contain personal patient information, it is the responsibility of the end-users to follow the local laws related to appropriate handling of personal data – for example HIPAA (USA) and GDPR (EU) – and to discard any information when required. Please refer to the relevant section in the Reference Guide regarding how personal data is stored and accessible within Simpleware Medical, to ensure your usage is compliant.

It is recommended to use Simpleware Medical within a hardware and/or network environment in which cyber security controls have been implemented including anti-virus and use of firewall. It is also recommended to ensure password protection for the workstation on which Simpleware Medical is running, and not to leave the workstation unattended.

Scripts written by the user and run through either the GUI or GUI-less, e.g. console mode, view in the software, should always be reviewed in the GUI version of the software before being used for any purpose. Scripts written for previous versions of Simpleware Medical should be re-validated by the user before being applied with the present UI-based software interface.

Users are also recommended to review the scripting API change log after upgrading to a newer version (available within the Scripting Help menu).

3.2 Accuracy

Simpleware Medical image processing and meshing algorithms are designed to use partial volume effects to improve surface accuracy. The reconstructed 3D surface typically has a maximal error of ½ of a voxel size.

Note: the accuracy of a model is dependent on the image resolution and the quality of the original scan. The accuracy of a model for simulation is also dependent on user requirements and choice of simulation software.

During surface reconstruction, error can be found near sharp edges, which are difficult to reconstruct when using any image-based meshing techniques. Excessive noise in scanned images can also affect surface reconstruction accuracy.

3.3 Other Ways of Viewing this Information and Emergencies

End-users of Simpleware Medical can request a free paper copy of this document at no charge. It will be delivered within 7 calendar days.

To do so, please contact simpleware@synopsys.com.

If information is needed in an emergency, please call +44(0)1392 642635.

In the event that you experience temporary unavailability of this document through the Synopsys website or to the Internet in general, or of your institutional access, we recommend temporarily suspending use of the software until access is restored, unless you have a paper copy of this document.

There are no foreseeable medical emergencies residual risks, adverse events, or side effects related to this device. If you believe that the device may have directly or indirectly contributed to a patient's injury or death, then please immediately contact simpleware@synopsys.com or call +44(0)1392 642635.

Please also contact the person in charge of communication with your Competent Authority, according to the procedures of your organization.

3.4 Dialog Symbols

Simpleware Medical uses a set of standard symbols (icons) when displaying information dialogs. The table below provides information about the severity of the risk associated with each type of symbol.

Symbol	Meaning
	<ul style="list-style-type: none"> Hint or Information. Risk level: low. The provided information may be useful but ignoring it has a very low risk of causing patient injury or death.
	<ul style="list-style-type: none"> Warning or potential usage error. Risk level: medium. The provided information is related to an operation that has succeeded but the user needs to be aware of details which may indicate undesirable side effects.
	<ul style="list-style-type: none"> Error/operation failure. Risk level: low. An operation has failed. Even though the outcome is negative, there is no risk of undesirable side effects.

3.5 Instructions for Use

3.5.1 Starting Simpleware Medical

The installation can be carried out by the user. There are no operational security options or critical dependencies on other applications required for installation of the software. Please see section **Warnings and Recommendations** of this document for more information on general cybersecurity requirements.

After installing the software on your PC, double click the Simpleware Medical V-2024.06 icon on your desktop. Alternatively, you can click on the Windows icon in the “taskbar” and navigate to **Synopsys Simpleware Medical > Simpleware Medical V-2024.06**.

The Reference Guide can be opened by clicking on the “Help” button on the Welcome page that appears after starting Simpleware Medical.

Simpleware Medical usage is controlled through a license key file which may either be node-locked or floating. Instructions for setting up both license options are described in the **Reference Guide**. All required software installers (product and licensing tools) and license keys can be downloaded from **SolvNetPlus** (<https://solvnet.synopsys.com>). Please note that only active licensees of Simpleware Medical will be able to get access to SolvNetPlus.

3.5.2 Supported Operating Systems

- Windows 10*
- Windows 11*†
- Windows Server 2016*
- Windows Server 2019*
- Windows Server 2022*
- Licensing tools for Linux 32-bit and 64-bit

*Only 64-bit versions of these operating systems are supported.

†Simpleware Medical is fully tested on this operating system.

3.5.3 Recommended System Requirements

Operating System	Hardware
<ul style="list-style-type: none"> Windows 11 Windows Server 2019 <p>Only 64-bit operating systems are supported.</p>	<ul style="list-style-type: none"> Processor: Intel Core i7 or equivalent Memory (RAM): 16 GB or higher** Graphics card: Workstation grade OpenGL 3.2 compatible graphics card with 1024 MB VRAM Screen resolution: 1920 x 1080; True color (32bit) or more Disk space: 100 GB <p>** As a guideline, allow twice as much RAM as your datasets size for volume rendering visualization and about 4 to 8 times for segmentation, surface mesh generation etc. The memory required depends heavily on the size and complexity of your image data and resultant model to be analyzed. Please contact Simpleware Support for further information with any specific queries.</p>

3.5.4 Uninstalling Simpleware Medical

Please refer to the steps for uninstalling programs on your particular version of Windows. In most cases, uninstalling can be completed by:

Using 'Add/remove programs' for the Windows settings app:

Step 1: In Windows, go to 'Start' (or press the "Windows" key).

Step 2: Search for 'Add and remove programs' and open the app.

Step 3: Search for 'Simpleware Medical'.

Step 4: Select the version and then select 'Uninstall'.

Using the Search in Windows taskbar:

Step 1: In Windows, go to 'Start' (or press the "Windows" key).

Step 2: Search for 'Simpleware Medical'.

Step 3: Under the App information on the right of the search results, expand the list using the down arrow.

Step 4: Select 'Uninstall'.

We also suggest deleting the software directory to ensure all data is removed, typically by going to Program Files and finding the folder for Simpleware Medical, if still present after uninstalling the software. Please also double-check that all data is removed from:

C:\Users\(or C:\Users\

C:\Windows\Temp folder (or where the user specified in the Simpleware Preferences)


4 Manufacturer Contact Details

Simpleware Medical V-2024.06

Manufactured in June 2024 by Synopsys Simpleware Product Group

Synopsys (Northern Europe) Ltd.
Bradninch Hall, Castle Street
Exeter, Devon, EX4 3PL
United Kingdom

Phone: +44(0)1392 642635
Email: simpleware@synopsys.com
SRN: GB-MF-000002376

2797

Simpleware Medical is a CE-marked product

UDI Number: (01)00860011545705(8012)V-2024.06

Synopsys International Ltd.
Floor 4, No. 2 Stemple Exchange, Blanchardstown Corporate Park, Dublin 15,
D15E4FN, Ireland
SRN: IE-AR-000001516

Synopsys Switzerland LLC
Thurgauerstrasse 40, 8050 Zürich, Switzerland
CHRN-AR-20003762

5 Current Availability as a CE-Marked Product

For more information on the current availability and plans for Simpleware Medical as a CE-marked product in your country, please contact simpleware@synopsys.com.

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Synopsys (Nordeuropa) Ltd.